

# Wiper/Chemical Compatability Chart

Class	Chemical	Nylon	Polyester	Polyester/Cellulose	Polypropylene	Polypropylene/Cellulose	Cotton	Rayon/Regenerated Cellulose
Acids	Acetic acid, 5%	R	R	L	R	L	L	R
	Acetic acid, 10%	L	R	L	R	L	L	R
	Acetic acid, Glacial	N	R	L	R	L	L	R
	Boric Acid	L	R	R	R	R	R	R
	Hydrochloric, 6N	N	L	L	R	R	R	N
	Hydrochloric, Conc.	N	N	N	L	N	N	N
	Hydrofluoric, 10%	N	R	L	L	L	N	N
	Hydrofluoric, 35%	N	R	N	N	N	N	N
	Nitric Acid, 6N	N	R	L	L	L	L	N
	Nitric Acid, Conc.	N	N	N	N	N	N	N
Alcohols	Sulfuric Acid, 6N	N	R	N	R	N	N	N
	Sulfuric Acid, Conc.	N	N	N	R	N	N	N
	Amyl Alcohol	R	R	L	R	L	L	R
	Benzyl Alcohol	L	R	N	R	N	N	R
	Butyl Alcohol	R	R	L	R	L	L	R
	Ethyl Alcohol, <80%	R	R	R	R	R	R	R
	Ethyl Alcohol, >80%	R	R	L	R	L	L	R
	Ethylene Glycol	R	R	L	R	L	L	R
	Glycerine (Glycerol)	R	R	R	R	R	R	R
	Isobutyl Alcohol	R	R	R	R	R	R	R
Bases	Isopropanol	R	R	L	R	L	L	R
	Methanol	L	R	T	R	T	T	R
	Methyl Cellosolve	R	R	L	R	L	L	L
	Propanol	R	R	R	R	R	R	R
Misc.	Ammonium Hydroxide, 6N	N	L	N	R	N	N	L
	Potassium Hydroxide, 6N	R	N	N	R	N	N	N
	Sodium Hydroxide, 6N	N	L	N	R	N	N	N
	Cottonseed Oil	R	T	R	R	R	R	T
	Hydrogen Peroxide (30%)	R	R	R	R	R	R	N
	Peanut Oil	R	R	R	R	R	R	T
	Petroleum Oils	T	R	R	R	R	R	R
	Sesame Oils	R	R	R	R	R	R	T
Misc.	Silicone Oils	R	R	R	R	R	R	R
	Turpentine	R	R	N	R	L	N	R

Class	Chemical	Nylon	Polyester	Polyester/Cellulose	Polypropylene	Polypropylene/Cellulose	Cotton	Rayon/Regenerated Cellulose
Solvents	Acetone	R	R	T	R	L	N	R
	Acetonitrile	R	R	N	R	L	N	R
	Amyl Acetate	R	R	N	R	L	N	R
	Aniline	R	R	N	R	L	N	R
	Benzene	R	R	N	R	L	N	R
	Bromoform	R	R	N	R	L	N	R
	Butyl Acetate	R	R	N	R	L	N	R
	Carbon Tetrachloride	R	R	N	R	L	N	R
	Cellosolve	R	R	N	R	L	N	L
	Chloroform	R	R	L	R	L	L	R
	Cyclohexane	R	R	L	L	L	L	R
	Cyclohexanone	R	R	N	L	L	N	R
	Diethyl Acetamide	R	R	T	T	T	T	R
	Dimethyl Formamide	R	R	N	R	L	N	L
	Dimethyl Sulfoxide (DMSO)	R	R	N	R	L	N	L
	Dioxane	R	R	N	R	L	N	L
	Ethyl Ether	R	R	T	T	T	T	R
	Ethylene Dichloride	R	R	N	L	L	N	R
	Formaldehyde	R	R	L	R	L	L	R
	Freon TF	R	R	R	R	R	R	R
	Gasoline	R	R	R	R	R	R	R
	Hexane	R	R	R	L	L	R	R
	Isopropyl Acetate	R	R	N	R	L	N	R
	Kerosene	R	R	R	R	R	R	R
	Methyl Acetate	R	R	R	R	R	N	R
	Methyl Ethyl Ketone (MEK)	R	R	N	R	L	N	R
	Methyl Isobutyl Ketone	R	T	N	R	L	N	R
	Methylene Chloride	L	R	L	R	L	L	R
	Nitrobenzene	R	R	N	N	L	N	L
	Pentane	R	R	R	R	R	R	R
	Perchloroethylene	R	T	N	L	R	N	R
	Pyridine	L	R	N	R	L	N	R
	Tetrahydrofuran	R	R	N	R	L	N	R
	Toluene	R	R	R	L	L	R	R
	Trichloroethane	R	T	L	R	L	L	R
	Trichloroethylene	R	R	R	L	L	R	R
	Triethylamine	R	R	N	L	N	N	R
	Xylene	R	R	N	L	L	N	R

Key:  
 R = Recommended  
 L = Limited Resistance  
 N = Not recommended  
 T = Test  
 (Testing before use is recommended)

1. Spectra/Por® Biotech Dialysis Membranes Cellulose Ester (CE) Regenerated Cellulose (RC). Spectrum Product Instruction Manual 420-10688-000Rev04. 17 January 2008 www.spectrapor.com/lit/420x10688x000.pdf  
 2. Chemical Compatibility. Spectrum Laboratories, Inc. 1995-2008.7 January 2008 www.spectrapor.com/cell/Compatibility.html  
 3. Filterware Chemical Resistance Database. NALGENE® Labware. 2008. 17 January 2008 www.nalgenelabware.com/techdata/Chemical/filter.asp

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at **Belgium:** fishersci.be **Denmark:** fishersci.dk  
 Germany: fishersci.de **Ireland:** fishersci.ie **Italy:** fishersci.it  
 Finland: fishersci.fi **France:** fishersci.fr **Netherlands:** fishersci.nl  
 Norway: fishersci.no **Portugal:** fishersci.pt **Spain:** fishersci.es  
 Sweden: fishersci.se **Switzerland:** fishersci.ch **UK:** fishersci.co.uk

