




**Thermo Scientific**  
Magnetic Bead Technology

# magnetic bead technology for better assay development

custom magnetic particles • immunoprecipitation • affinity purification • ChIP

**Thermo**  
SCIENTIFIC


exceptional bead technology  
for **protein biology**




Streptavidin

**Thermo Scientific™ Magnetic Beads** are available for


immunoprecipitation and affinity purification as well as custom magnetic particle creation. The high-performance, iron oxide, superparamagnetic particles are validated and optimized for use with high-throughput magnetic platforms, such as the Thermo Scientific™ KingFisher™ Duo and Flex Instruments. Samples can be analyzed by Western blotting or on Thermo Scientific™ Mass Spectrometers for quantitation of low abundant targets.



Identify and  
quantitate



Identify and quantitate  
your analyte using  
Thermo Scientific Mass  
Spectrometers and  
labeling reagents



detect

Use our highly sensitive  
substrates for Western  
blot detection

lyse → target → isolate →



Prepare cells, tissue or  
fluid samples with Thermo  
Scientific™ Lysis Reagents




Isolate your target with one  
of our magnetic beads:

- **Activated magnetic beads** for creating custom affinity resins .. pg. 2
- **Magnetic Protein A/G/L** for immunoprecipitation ..... pg. 4
- **Validated kits** for immunoprecipitation, ChIP or RNA-pull downs ..... pg. 16
- **Resins** for protein enrichment ..... pg. 20



Choose to run a few  
samples on the bench or  
create a high-throughput  
assay with KingFisher  
Instruments




# immobilize your ligand

for custom magnetic assays

## Magnetic Protein Immobilization Beads


**Thermo Scientific™ Pierce™ NHS-Activated Magnetic Beads** enable covalent, amine-based conjugation of proteins to magnetic beads in a simple mix-and-go format for use in custom affinity purification experiments. The activated magnetic beads contain *N*-hydroxysuccinimide (NHS) functional groups that react with primary amines, forming stable amide linkages. Once they are covalently attached, the immobilized proteins are highly resistant to leaching from the bead surface. When prepared beads are used in experiments, nonspecific binding is negligible because nonreacted NHS-ester groups are thoroughly blocked during the coupling procedure.







**Figure 1. Significantly better coupling capacity with Thermo Scientific Pierce NHS-Activated Magnetic Beads.** Rabbit IgG (1 mg/mL) was coupled in PBS for two hours at pH7.2 to 3mg each of Pierce NHS-Activated Magnetic Beads and NHS Mag Sepharose™ Beads (GE Life Sciences). Negative control beads (-) were prepared by quenching or blocking using respective manufacturer protocols. Bound protein was measured using the Thermo Scientific™ Pierce™ 660nm Protein Assay by subtracting the amount of protein in the flow-through from the amount loaded. The Pierce NHS-Activated Magnetic Beads coupled more than four times as much protein as the equivalent amount of NHS Mag Sepharose Beads.

### Highlights

- **High capacity** – at least four times greater binding capacity than NHS-activated magnetic beads from other suppliers
- **Easy to use** – immobilize in a simple one-step reaction with minimal hands-on time
- **Safe** – no hazardous chemicals (e.g., sodium cyanoborohydride and cyanogen bromide) needed
- **Ligand compatible** – use with nearly any primary amine-containing compound or affinity ligand to immobilize
- **Low nonspecific binding** – the bead surface is pre-blocked and any nonreacted NHS-ester groups are fully quenched
- **Protocol compatible** – protein coupling to the beads and downstream applications can be performed manually or by automation (e.g., KingFisher Instruments)



**Figure 2. Reaction scheme for conjugation of protein onto Thermo Scientific Pierce NHS-Activated Magnetic Beads.**



**Table 1. Properties of Thermo Scientific Pierce NHS-Activated Magnetic Beads.**

<b>Composition</b>	<i>N</i> -hydroxysuccinimide (NHS) functional groups on a blocked magnetic bead surface
<b>Magnetization</b>	Superparamagnetic (no magnetic memory)
<b>Mean Diameter</b>	1 $\mu$ m (nominal)
<b>Density</b>	2.0g/cm <sup>3</sup>
<b>Bead Concentration</b>	10mg/mL in DMAC
<b>Binding Capacity</b>	$\geq 26\mu$ g of rabbit IgG/mg of beads



### Ordering Information

Product #	Description	Pkg. Size
<b>88826</b>	<b>Pierce NHS-Activated Magnetic Beads</b> Sufficient for: Binding $\geq 26\mu$ g of rabbit IgG/mg of beads	1 mL
<b>88827</b>	<b>Pierce NHS-Activated Magnetic Beads</b> Sufficient for: Binding $\geq 26\mu$ g of rabbit IgG/mg of beads	5 mL

**Figure 3. Effective cell separation.** Pierce NHS-Activated Beads were coated with Thermo Scientific™ Stage-specific Embryonic Antigen 1 (SSEA-1) mouse IgM or rabbit IgG as a negative control. The beads were incubated with a 50:50 co-culture of F9 mouse embryonal carcinoma cells (SSEA-1 positive) and NIH 3T3 cells (SSEA-1 negative) for 20 minutes at 4°C. The beads were collected on a magnetic stand and the unbound cell fraction was evaluated by flow cytometry using anti-SSEA-1 mouse IgM and goat anti-mouse IgM-fluorescein. **Panel A** shows F9 cells were selectively depleted with anti-SSEA-1-coated NHS-activated magnetic beads. **Panel B** shows neither cell type bound to the negative control (rabbit IgG-coated NHS magnetic beads). Both the bead-bound and unbound cell fractions were cultured for 24 hours, fixed and then stained with mouse anti-SSEA-1 antibody, Thermo Scientific™ DyLight™ 488 conjugated goat anti-mouse IgM and Hoechst™ nuclear stain. Cells were visualized using the Thermo Scientific™ ToxInsight™ Platform. **Panel C** shows that the SSEA-1 antibody coated magnetic NHS beads effectively separated the F9 cells from the NIH 3T3 cells. As expected, the rabbit-IgG-coated magnetic NHS beads (**Panel D**) did not bind either cell type, and the corresponding unbound fraction contained a 50:50 ratio of both cell types.

# build your own immunoprecipitation assay

## Magnetic Immunoprecipitation Beads

The Thermo Scientific™ Pierce™ Protein A/G Magnetic Bead is a single particle that is compatible with all commonly used antibodies for immunoprecipitation (IP). These beads are coated with genetically engineered Protein A/G, a recombinant protein that combines the IgG binding domains of both Protein A and Protein G. This combination enables the capture of antibodies from a wider range of species and isotypes than either protein alone. Using our crosslinker chemistry, you can immobilize an antibody onto the magnetic particle and prevent IgG contamination in your immunoprecipitated sample.




Figure 4. Diagram of Thermo Scientific Pierce Protein A/G Magnetic Beads. The magnetic particles are 1 µm in diameter and are specially manufactured with two layers of magnetite and encapsulation. Recombinant Protein A/G is coupled to the bead surface.

### Highlights

- **Compatible** – one magnetic bead type that can capture most primary antibodies
- **Fast** – immunoprecipitating in as few as 30 minutes helps reduce nonspecific binding and improves the capture of transient protein complexes
- **Clean** – immobilize your antibody to prevent contamination in your eluate
- **Resistant** – no leaching of Protein A/G in the presence of detergents, low pH buffers or common mass spectrometry solvents
- **Efficient** – immunoprecipitate with half the recommended volume of magnetic particles compared to other magnetic beads

Table 2. Properties of Thermo Scientific Pierce Magnetic Protein A/G Beads.

Composition	Magnetite-coated polymeric beads blocked and covalently coated with a monolayer of recombinant Protein A/G
Mean Diameter	1 µm (nominal)
Density	2.0g/cm <sup>3</sup>
Bead Concentration	10mg/mL in water with sodium azide
Binding Capacity	55 to 85µg of rabbit IgG/mg magnetic beads








Figure 5. Thermo Scientific Pierce Protein A/G Magnetic Beads isolate significantly more IgG from rabbit and mouse serum with less background than other brands of Protein A and Protein G magnetic particles. Using a KingFisher Flex Instrument with a 96 deep well plate, IgG was purified from 5mg of rabbit and mouse serum using 50µL of Pierce Protein A/G Magnetic Beads, Dynabeads™ Protein A or G (Life Technologies), or PureProteome™ Protein A or G Beads (Millipore). The beads were washed with Tris-buffered saline containing 0.05% Tween™ -20 (TBST), incubated for one hour with serum diluted in TBST, washed three times, and then eluted with 0.1M glycine, pH 2.8 for 10 minutes at room temperature. The eluates were resolved by SDS-PAGE and stained with Thermo Scientific™ Imperial™ Protein Stain. **Panel A:** Rabbit serum; **Panel B:** Mouse serum. The IgG heavy chain bands were quantified by densitometry. The values for each set of duplicate bands were averaged and expressed as a percentage of the average for the Pierce Protein A/G Magnetic Beads.



**Figure 6.** The rabbit IgG binding capacity of Thermo Scientific Pierce Protein A/G Magnetic Beads is approximately four times greater than that of other Protein A and Protein G beads. Pierce Protein A/G Magnetic Beads or Dynabeads Protein A or Protein G (Life Technologies) were added to a 96 deep-well plate (1 mg beads per well). Using the Thermo Scientific™ KingFisher™ 96 Instrument, the beads were incubated for one hour with varying amounts of purified rabbit IgG (20 to 200µg). After binding, the samples were eluted at 96°C with SDS-PAGE reducing sample buffer. Binding was calculated using the Thermo Scientific™ Pierce™ BCA Protein Assay.



**Figure 7.** The Thermo Scientific Pierce Protein A/G Magnetic Beads effectively immunoprecipitate cell cycle proteins Cyclin D, Cyclin E, Cyclin B and Cdk1. U2OS (human osteosarcoma) cells were synchronized at G<sub>0</sub> followed by growth in 20% fetal bovine serum for 4, 6 and 18 hours before harvest. The cells were lysed in IP lysis/wash buffer, and 0.75mg of lysate (per sample) was incubated with anti-Cyclin D (rabbit polyclonal), anti-Cyclin E (mouse IgG<sub>1</sub>), anti-Cyclin B (mouse IgG<sub>1</sub>) or anti-Cdk1 (rabbit polyclonal) antibodies overnight at 4°C. The Pierce Protein A/G Magnetic Beads were added (50µL each) to a 96 deep-well plate and immunoprecipitations were performed using the KingFisher Flex Instrument. Eluted sample volumes of 5µL, 10µL and 20µL were resolved by SDS-PAGE and analyzed by Western blot.



### Ordering Information

Product #	Description	Pkg. Size
88802	<b>Pierce Protein A/G Magnetic Beads</b> Sufficient for: Binding 55 to 85µg rabbit IgG/mg beads.	1 mL
88803	<b>Pierce Protein A/G Magnetic Beads</b> Sufficient for: Binding 55 to 85µg rabbit IgG/mg beads.	5 mL

# build your own immunoprecipitation assay

## Magnetic Immunoprecipitation Beads


**Thermo Scientific™ Pierce™ Protein A Magnetic Beads** are used for immunoprecipitating antigens from cell or tissue extracts as well as purifying antibody from serum, cell culture supernatant or ascites fluid. Protein A can bind to antibodies from many different species, including mouse, human, rabbit, pig, dog and cat. The protocol for Pierce Protein A Beads is optimized for high recovery and high purity of isolated antibodies or antigens.

### Highlights


- **Low nonspecific binding** – stable, pre-blocked beads provide clean purification product
- **Consistency** – magnetic beads eliminate resin loss and provide for more efficient separation of solutions than traditional IP methods that use only microcentrifuge tubes
- **Compatibility** – beads are compatible with manual and automated applications (e.g., KingFisher Instruments)

Table 3. Properties of Thermo Scientific Pierce Magnetic Protein A Beads.

Composition	Magnetite-coated polymeric beads blocked and covalently coated with a monolayer of recombinant Protein A
Mean Diameter	1 μm (nominal)
Density	2.0g/cm <sup>3</sup>
Bead Concentration	10mg/mL in water with sodium azide
Binding Capacity	≥40μg of rabbit IgG/mg of beads; ≥400μg of rabbit IgG/mL of beads




**Figure 8.** The human and rabbit IgG binding capacities of Thermo Scientific Pierce Protein A Magnetic Beads are approximately 2-fold higher than LifeTech Dynabead Protein A. Pierce Protein A Magnetic Beads or Dynabeads A (Life Technologies) were added to a 96 deep-well plate (100μL beads per well). Using the KingFisher Flex Instrument, the beads were incubated for one hour with 400μg purified human or rabbit IgG. Binding was calculated using the Pierce BCA Protein Assay by subtracting the amount of IgG in the flow-throughs from the IgG loaded.




**Figure 9.** Thermo Scientific Pierce Protein A Magnetic Beads exhibit low nonspecific binding. Using a KingFisher Flex Instrument with a 96 deep-well plate, IgG was purified from 2mg of rabbit serum using 50μL of Pierce Protein A Magnetic Beads. The beads were incubated for one hour with serum diluted in phosphate-buffered saline containing 0.025% Tween-20 (PBST), washed twice with PBST and once with water, and then eluted with 0.1M glycine, pH 2.0 for 10 minutes at room temperature. The eluates were resolved and stained with Imperial Protein Stain. No serum proteins other than antibody heavy and light chains were detected in the eluted sample.






A.



B.



### Ordering Information

Product #	Description	Pkg. Size
88845	Pierce Protein A Magnetic Beads	1 mL
88846	Pierce Protein A Magnetic Beads	5 mL

**Figure 10. Better immunoprecipitation results with Thermo Scientific Pierce Protein A Magnetic Beads.** PP2A antibody (5 $\mu$ g) was incubated overnight at 4 $^{\circ}$ C with 0.5mg of A549 cell lysate. Using the KingFisher Flex Instrument, 50 $\mu$ L each of Pierce Protein A Magnetic Beads, Dynabeads A (Life Technologies), Protein A Magnetic Sepharose (GE Life Sciences) and PureProteome Protein A Beads (EMD/Millipore) were added to 96 deep-well plates. The beads were incubated for one hour with the antigen/antibody complex at room temperature, washed twice in phosphate-buffered saline containing 0.05% Tween-20, washed once in water and then eluted in 0.1M glycine, pH 2.0. Samples were resolved by SDS-PAGE and analyzed for Western blot for PP2A (**Panel A**) and by silver stain for nonspecific binding (**Panel B**). The Pierce Protein A Magnetic Beads were found to have higher yield of PP2A than other Protein A beads. Nonspecific binding was negligible for all beads tested.

# build your own immunoprecipitation assay

## Magnetic Immunoprecipitation Beads

### Thermo Scientific™ Pierce™ Protein G Magnetic Beads

are high-capacity and high-throughput affinity particles for antibody purification and immunoprecipitation methods using manual or robotic magnetic separators. Protein G can bind to antibodies from many different species, including mouse, human, rabbit, cow, goat and sheep.

#### Highlights

- **High IP efficiency** – highest antibody yield
- **Low nonspecific binding** – stable, pre-blocked beads provide clean purification product
- **Assay consistency** – magnetic beads eliminate resin loss
- **High throughput** – compatible with manual and automated applications (e.g., KingFisher Instruments)

Table 4. Properties of Thermo Scientific Pierce Magnetic Protein G Beads.

Composition	Magnetite-coated polymeric beads blocked and covalently coated with a monolayer of recombinant Protein G
Mean Diameter	1µm (nominal)
Density	2.0g/cm <sup>3</sup>
Bead Concentration	10mg/mL in water with sodium azide
Binding Capacity	≥60µg of rabbit IgG/mg of beads; ≥600µg of rabbit IgG/mL of beads




Figure 11. The human and rabbit IgG binding capacities of Thermo Scientific Pierce Protein G Magnetic Beads are approximately 3-fold higher than Dynabeads Magnetic Beads. Pierce Protein G Magnetic Beads and Dynabeads G (Life Technologies) were added to a 96 deep-well plate (100µL beads per well). Using the KingFisher Flex Instrument, the beads were incubated for one hour with 400µg purified human or rabbit IgG. Binding was calculated using the Pierce BCA Protein Assay by subtracting the amount of IgG in the flow-throughs from the IgG loaded.






Figure 12. Thermo Scientific Pierce Protein G Magnetic Beads exhibit low nonspecific binding. Using a KingFisher Flex Instrument with a 96 deep-well plate, IgG was purified from 2mg of rabbit serum using 50µL of Pierce Protein G Magnetic Beads. The beads were incubated one hour with serum diluted in phosphate-buffered saline containing 0.025% Tween-20 (PBST), washed twice with PBST and once with water, and then eluted with 0.1M glycine, pH 2.0 for 10 minutes at room temperature. The eluates were resolved and stained with Imperial Protein Stain. No serum proteins other than antibody heavy and light chains were detected in the eluted sample.




Panel A.



### Ordering Information

Product #	Description	Pkg. Size
88847	Pierce Protein G Magnetic Beads	1mL
88848	Pierce Protein G Magnetic Beads	5mL

Panel B.



**Figure 13. Better immunoprecipitation results with Thermo Scientific Pierce Protein G Magnetic Beads.** EGF Receptor antibody (5µg) was incubated overnight at 4°C with 0.75mg of A431 cell lysate. Using the KingFisher Flex Instrument, 25µL each of Pierce Protein G Magnetic Beads, Dynabeads G (Life Technologies), Protein G Magnetic Sepharose (GE Life Sciences) and PureProteome Protein G Beads (EMD/Millipore) were added to 96 deep-well plates. The beads were incubated for one hour with the antigen/antibody complex at room temperature, washed twice in phosphate-buffered saline containing 0.05% Tween-20, washed once in water and then eluted in 0.1M glycine, pH 2.0. Samples were resolved by SDS-PAGE and analyzed by Western blot for EGFR (**Panel A**) and by silver stain for nonspecific binding (**Panel B**). The Pierce Protein G Magnetic Beads were found to have higher yield of EGFR than other Protein G beads. Nonspecific binding was negligible for all beads tested.


# build your own immunoprecipitation assay

## Magnetic Immunoprecipitation Beads

**Thermo Scientific™ Pierce™ Protein L Magnetic Beads** are ideal for selective isolation of antibodies possessing kappa light chains. Protein L selectively binds mouse and human antibodies through kappa light chains and is commonly used to purify monoclonal antibodies in cell culture supernatants supplemented with bovine serum as Protein L does not bind bovine IgG. Protein L can bind a broader range of Ig classes than Protein A or Protein G, including IgG, IgM, IgA, IgE and IgD. Protein L binds strongly to human (kappa I, III and IV only), mouse (kappa I only), rat and pig immunoglobulins. It binds weakly to rabbit immunoglobulins and does not bind to immunoglobulins from bovine, goat or sheep. Single-chain variable fragments (scFv) and Fab fragments also bind to Protein L.

### Highlights



- **Selective** – ideal for selective purification of human and mouse antibodies that have kappa light chains
- **Low nonspecific binding** – stable, pre-blocked beads provide clean purification of antibody
- **Compatibility** – beads are compatible with manual and automated applications (e.g., KingFisher Instruments)



**Figure 14. Human and mouse IgG binding capacity curves for Thermo Scientific Pierce Protein L Magnetic Beads.** Pierce Protein L Magnetic Beads were added to a 96 deep-well plate (100µL beads per well). Using the KingFisher Flex Instrument, the beads were incubated for one hour with purified human or mouse IgG (amounts shown in graphs). Binding was calculated using the Pierce BCA Protein Assay by subtracting the amount of IgG in the flow-throughs from the IgG loaded.

**Table 5. Properties of Thermo Scientific Pierce Magnetic Protein L Beads.**

<b>Composition</b>	Magnetite-coated polymeric beads blocked and covalently coated with a monolayer of recombinant Protein L
<b>Mean Diameter</b>	1µm (nominal)
<b>Density</b>	2.0g/cm <sup>3</sup>
<b>Bead Concentration</b>	10mg/mL in water with sodium azide
<b>Binding Capacity</b>	≥110µg of human IgG/mg of beads; ≥1.1mg of human IgG/mL of beads



### Ordering Information

Product #	Description	Pkg. Size
88849	Pierce Protein L Magnetic Beads	1 mL
88850	Pierce Protein L Magnetic Beads	5 mL

**Figure 15. Thermo Scientific Pierce Protein L Magnetic Beads isolate more mouse IgG from cell culture supernatant than other suppliers' Protein L magnetic beads.** Using a KingFisher Flex Instrument with 96 deep-well plates, IgG was purified from cell culture supernatant using 50µL each of Pierce Protein L Magnetic Beads, and Protein L Magnetic Beads from BioVision and Bioclone. The beads were incubated for one hour with undiluted cell culture supernatant containing 0.025% Tween-20, washed twice with PBST and once with water, and then eluted with 0.1M glycine, pH 2.0 for 10 minutes at room temperature. The eluates were resolved and stained with the Thermo Scientific™ Pierce™ Silver Stain Kit. Note that binding of mouse IgG with Protein L only occurs when kappa light chains are present. All of the beads were found to have negligible nonspecific binding.

# select an easy-to-use validated kit

## Magnetic Immunoprecipitation (IP) Kits

**Thermo Scientific™ Pierce™ Magnetic IP/Co-IP Kits** are optimized to isolate protein complexes from biological samples. Each kit contains all the required buffers and beads validated to deliver the best results. Three versions of the kit are available to perform a classic IP, crosslink IP or direct IP.



### Kit Highlights

- Compatible with any antibody
- Faster immunoprecipitations (IPs) for less background
- Easily capture transient protein complexes
- No antibody contamination in your eluted sample
- Simple handling with no sample loss
- Validated for automated protocols using KingFisher Instruments

**The Thermo Scientific™ Pierce™ Classic Magnetic IP/Co-IP Kit** uses high binding capacity Pierce Magnetic Protein A/G Beads to deliver clean and consistent co-immunoprecipitations (co-IP) with any common antibody. Antibodies are not linked to the resin and will co-elute with your antigen.

### Select this version:


- For the highest antigen yield
- If antibody contamination is not a concern



**Figure 16. The Thermo Scientific Pierce Classic Magnetic IP/Co-IP Kit immunoprecipitates Grp94 with higher yield than other Protein A and Protein G beads.** MOPC (mouse myeloma) cells were lysed in RIPA buffer and 0.75mg of lysate was incubated with Grp94 antibody (rat IgG<sub>2a</sub>) overnight at 4°C. Using the KingFisher Flex Instrument, 50µL of Pierce Protein A/G Magnetic Beads and 50µL each of Protein A and Protein G beads from Life Technologies, EMD/ Millipore and GE Life Sciences were added to a 96 deep-well plate. The eluates were resolved by SDS-PAGE and analyzed by Western blot for Grp94.

### Ordering Information

Product #	Description	Pkg. Size
88804	<b>Pierce Classic Magnetic IP/Co-IP Kit</b> <i>Sufficient for: 40 IP reactions using 25µL of beads</i> Contains: Pierce Protein A/G Magnetic Beads, 1mL Pierce IP Lysis/Wash Buffer, 2 x 50mL Lane Marker Sample Buffer (5X), 5mL Elution Buffer, 5mL Neutralization Buffer, 0.5mL	40-rxn kit







**The Thermo Scientific™ Pierce™ Crosslink Magnetic IP/Co-IP Kit** uses crosslinkers to immobilize your primary antibody to Protein A/G. This prevents antibody contamination in your eluted sample and eliminates antibody interference in Western blot and mass spec applications.

**Select this version:**

- To eliminate antibody contamination that interferes with downstream detection
- To properly orient your antibody





**Figure 17. The Thermo Scientific Pierce Crosslink Magnetic IP/Co-IP Kit immunoprecipitates PP2A without antibody contamination and with negligible background.** PP2A antibody (5µg) was coupled to Pierce Protein A/G Magnetic Beads without DSS crosslinking (traditional IP) and with DSS crosslinking (crosslink IP). The beads were incubated with 0.5mg of A549 cell lysate for one hour at room temperature on the KingFisher Flex Instrument. PP2A was eluted from the beads with elution buffer for five minutes at room temperature and then neutralized with neutralization buffer. The eluates, antibody control (Ab) and flow-through (FT) were resolved by SDS-PAGE and analyzed by Western blot for PP2A. The antibody-crosslinked Pierce Protein A/G Magnetic Beads effectively immunoprecipitated PP2A without antibody contamination whereas the traditional IP method resulted in significant antibody contamination in the eluate.

**The Thermo Scientific™ Pierce™ Direct Magnetic IP/Co-IP Kit** uses Pierce NHS-Activated Magnetic Beads to immobilize your primary antibody directly to the bead surface. This method is independent of antibody species and prevents antibody contamination in your eluted sample.

**Select this version:**

- For non-traditional antibodies that do not bind Protein A or Protein G
- To eliminate antibody contamination



**Figure 18. Better immunoprecipitation results with Thermo Scientific Pierce NHS-Activated Magnetic Beads.** Anti-FOXP2 antibody (5µg) was coupled to 25µL of Pierce NHS-Activated Magnetic Beads and an equivalent amount of NHS Mag Sepharose (GE Life Sciences). The two sets of prepared beads were then used to immunoprecipitate FOXP2 from 0.5mg aliquots of the same 293T (human epithelial kidney) cell lysate. The eluates were resolved by SDS-PAGE and analyzed by Western blot for FOXP2.

**Ordering Information**

Product #	Description	Pkg. Size
<b>88805</b>	<b>Pierce Crosslink Magnetic IP/Co-IP Kit</b> Sufficient for: 40 IP reactions using 25µL of beads Contains: Pierce Protein A/G Magnetic Beads, 1mL IP Lysis/Wash Buffer, 2 x 50mL Coupling Buffer (20X), 25mL DSS Crosslinker, 8 x 2mg Lane Marker Sample Buffer (5X), 5mL Elution Buffer, 10mL Neutralization Buffer, 1mL Lane Marker Sample Buffer (5X), 5mL	40-rxn kit

**Ordering Information**

Product #	Description	Pkg. Size
<b>88828</b>	<b>Pierce Direct Magnetic IP/Co-IP Kit</b> Sufficient for: 40 IP reactions using 25µL of beads Contains: Pierce NHS-Activated Magnetic Beads, 1mL IP Lysis/Wash Buffer, 2 x 50mL Elution Buffer, pH 2.0, 5mL Lane Marker Sample Buffer, Non-reducing, (5X), 5mL Neutralization Buffer, pH 8.5, 0.5mL 0.67M Borate Buffer, 1mL BupH™ Borate Buffer Pack, 1 pack Quenching Buffer, 25mL	40-rxn kit

# select an easy-to-use validated kit

## Magnetic Anti-HA IP Kit

The Thermo Scientific™ Pierce™ HA-Tag Magnetic IP/Co-IP Kit provides a simple and fast method to study protein interactions. The high affinity anti-HA antibody-coupled magnetic beads enables immunoprecipitation (IP) of HA-tagged proteins or co-immunoprecipitation (co-IP) of their interacting partners without antibody contamination.

### Highlights

- **Specific** – immunoprecipitate only HA-tagged proteins and their interactors
- **Validated** – IP/co-IP kit includes a positive control lysate
- **Robust** – compatible with common tissue culture cell lysates
- **Convenient and easy** – complete kit includes all necessary reagents to perform 40 reactions

Table 6. Properties of Thermo Scientific Pierce Anti-HA Magnetic Beads.

Composition	Magnetite-coated polymeric beads blocked and covalently coated with a mouse monoclonal IgG1 anti-HA antibody
Mean Diameter	1µm (nominal)
Density	2.0g/cm <sup>3</sup>
Bead Concentration	10mg/mL in water with sodium azide
Binding Capacity	> 10µg of HA-tagged protein/mg beads (70kDA)




Figure 19. Significantly higher binding capacity with Thermo Scientific Pierce Anti-HA Magnetic Beads. HA-tagged protein (HA-ERK-GST, 400µg in PBS) was incubated with 100µL each of Thermo Scientific™ Pierce™ Anti-HA Magnetic Beads or MBL Anti-HA Magnetic Beads for one hour at room temperature. Bound HA-tagged protein was measured using the Pierce BCA Assay. Pierce Anti-HA Magnetic Beads pulled down more than four times as much protein as the equivalent amount of MBL Anti-HA Magnetic Beads.







Figure 20. Better Immunoprecipitation results from *E. Coli* lysates. Using a KingFisher Flex Instrument with 96 deep-well plates, 25µL each of Pierce Anti-HA Magnetic Beads, Anti-HA-tag Magnetic Beads (MBL International Corp.) and SPHERO™ Rabbit Anti-HA Magnetic Beads (Spherotech Inc.) were used to immunoprecipitate GST-PI3K-SH2-HA from 50µg of *E. coli* lysate in duplicate. Captured protein was eluted with 0.1M glycine, pH 2.0, and then resolved by SDS-PAGE and analyzed by Western blot for the HA-tagged protein.







**Figure 21. Better immunoprecipitation of protein expressed *in vitro*.** Using the KingFisher Flex Instrument, 25 $\mu$ L each of Pierce Anti-HA Magnetic Beads and Anti-HA-tag Magnetic Beads (MBL International Corp.) were added to 96 deep-well plates. The beads were incubated for one hour with HA-tagged BAD protein, expressed using the Thermo Scientific™ 1-Step High-Yield *in vitro* Translation Kit. After one hour at room temperature, samples were washed twice in phosphate-buffered saline containing 0.05% Tween-20, washed once in water and then eluted in 0.1M glycine, pH 2.0. Samples were resolved by SDS-PAGE and analyzed by Western blot for HA (**Panel A**) and by silver stain for nonspecific binding (**Panel B**). HA-BAD fusion protein yield was highest with Pierce Anti-HA Magnetic Beads compared to other anti-HA beads. Nonspecific binding was negligible and there was co-immunoprecipitation of the protein 14-3-3.



**Figure 22. Better co-IP results with Thermo Scientific Pierce Anti-HA Magnetic Beads.** Serine phosphorylation of BAD is associated with 14-3-3 binding and inhibition of BAD-induced cell death. Using a magnetic stand, 50 $\mu$ L each of Pierce Anti-HA Magnetic Beads and Anti-HA-tag Magnetic Beads (MBL International Corp.) were added to microcentrifuge tubes. The beads were incubated for one hour at room temperature with HA-tagged BAD expressed in the 1-Step Human High-Yield IVT Kit. After incubation, the beads were washed twice in phosphate-buffered saline containing 0.05% Tween-20, washed once in water and then eluted in 30% acetonitrile/0.5% formic acid. Samples were dried down in a speedvac and brought back up in 50 $\mu$ L of reducing SDS-PAGE samples buffer. One half of the reconstituted eluate was resolved by SDS-PAGE and analyzed by Western blot for 14-3-3. Pierce Anti-HA Magnetic Beads were found to have higher yield of 14-3-3 than the other anti-HA beads. Nonspecific binding was negligible.

### Ordering Information

Product #	Description	Pkg. Size
88836	Pierce Anti-HA Magnetic Beads	1 mL
88837	Pierce Anti-HA Magnetic Beads	5 mL
88838	Pierce HA-Tag Magnetic IP/Co-IP Kit Sufficient For: 40 IP reactions using 25 $\mu$ L of anti-HA magnetic beads Kit Contents: HA-tagged Positive Control (Product # 26180X), 500 $\mu$ L Application Set (Product # 88838X): Pierce Anti-HA Magnetic Beads, 0.65mL Pierce IP Lysis/Wash Buffer, 2 $\times$ 50mL, pH 7.4 Lane Marker Sample Buffer, Non-reducing, (5X), 5mL, pH 6.8 Elution Buffer, 5mL, pH 2.0 Neutralization Buffer, 1mL, pH 8.5	40-rxn kit

# select an easy-to-use validated kit

## Magnetic ChIP Kit

The Thermo Scientific™ Pierce™ Magnetic ChIP Kit provides a simple, fast and reproducible method to perform chromatin immunoprecipitation (ChIP) assays to capture a snapshot of specific protein-DNA interactions as they occur in living cells and then quantitate the interactions using PCR.

The Pierce Magnetic ChIP Kit contains sufficient reagents to perform 30 ChIP assays with appropriate controls using an optimized protocol. The blocked Pierce Protein A/G Magnetic beads used in this kit provide high binding capacity, low nonspecific background and flexibility of antibody species. These beads can be used manually with a magnetic stand as well as with automated platforms such as KingFisher Instruments. This kit provides reagents and a method to capture protein-DNA interactions *in vivo* allowing relative protein binding events to be monitored under different conditions and/or treatments. ChIP-validated and quality-guaranteed antibodies are also available for use with the Pierce Magnetic ChIP Kit.

### Highlights

- **Simple and fast** – obtain purified DNA ready for PCR in about eight hours
- **Efficient and reproducible** – micrococcal nuclease digestion and nuclear lysis are highly optimized
- **Sensitive** – obtain results with as little as  $1 \times 10^4$  cells
- **Low nonspecific background** – Pierce Protein A/G Magnetic beads are blocked in a non-DNA-containing reagent to minimize background
- **Complete** – optimized positive control reagents are included: RNA polymerase II antibody and GAPDH promoter PCR primers

### ChIP Assay Procedure

Total time: 8 Hours









Figure 23. Overview of the Thermo Scientific Pierce Magnetic ChIP Kit protocol. The Pierce Magnetic ChIP assay protocol and reagents provide an optimized system for performing chromatin crosslinking, cell lysis, IP and target protein recovery about eight hours.



**Figure 24.** The Thermo Scientific Pierce Magnetic ChIP Kit has a broad range of sensitivity. A431 lung carcinoma cells were crosslinked using a final concentration of 1% formaldehyde for 10 minutes. ChIP assays were performed with the Pierce Magnetic ChIP Kit to determine binding of RNA polymerase II to the proximal GAPDH promoter. Quantitative real-time PCR data was obtained with a Bio-Rad iQ5™ Thermocycler. Each column represents the fold enrichment of the RNA polymerase II over the normal rabbit IgG using the noted starting cell number (i.e., chromatin from that number of cells).



**Figure 25.** Greater fold enrichment than other kits. LNCaP prostate carcinoma cells were cultured in RPMI-1640 containing 10% FBS for 24 hours. Half of the cultures plated were treated for 16 hours with 5μM camptothecin, a drug that inhibits DNA topoisomerase I. Crosslinking was achieved using a final concentration of 1% formaldehyde in the media for 10 minutes. ChIP assays were performed according to the manufacturers' protocols to determine binding of p53 to a 1.5-kb region of the CDKN1A (p21) promoter. Quantitative real-time PCR data was obtained with a Bio-Rad iQ5 Thermocycler. The Pierce Magnetic ChIP Kit has been optimized to isolate even large DNA fragments.



**Figure 26.** Profiling multiple transcription factors binding to the MYC promoter using the Thermo Scientific Pierce Magnetic ChIP Kit. A431 lung carcinoma cells were cultured in DMEM containing 10% FBS for 24 hours. Following a 24-hour serum withdrawal, half of the cultures plated were treated with 100ng/mL EGF for 10 minutes. Crosslinking was achieved using a final concentration of 1% formaldehyde in the media for 10 minutes. ChIP assays were performed with the Pierce Magnetic ChIP Kit to determine binding of phosphorylated-STAT3, CBP, acetylated-CBP, histone H3, and RNA polymerase II to the proximal MYC promoter. Primary antibody amounts were determined empirically. Quantitative real-time PCR data was obtained with a Bio-Rad iQ5 Thermocycler. Graph represents the signal relative to the total input of chromatin.




### Ordering Information

Product #	Description	Pkg. Size
26157	<b>Pierce Magnetic ChIP Kit</b> Sufficient reagents to perform 30 ChIP Reactions	30-rxn kit
26162	<b>ChIP-grade Protein A/G Magnetic Beads</b> Formulation: Magnetite- and protein-coated polymer beads at 10mg/mL in water with sodium azide. Sufficient For: Use in approx. 250 typical ChIP assays	5mL

# select an easy-to-use validated kit

## Magnetic RNA-Protein Pull-Down Kit


**The Thermo Scientific™ Pierce™ Magnetic RNA-Protein Pull-Down Kit** provides researchers with a streamlined, robust method for enrichment and identification of RNA binding proteins. This method uses RNA probes labeled at the 3'-end with desthiobiotin and magnetic streptavidin particles. The complete kit contains sufficient reagents for 20 RNA-labeling reactions and 20 RNA-protein pull-down assays. Both synthetic RNA or *in vitro* transcribed RNA can be labeled with desthiobiotin. RNA-binding proteins are then enriched from cellular/tissue lysates or from *in vitro* translated protein preps. RNA-binding proteins are detected using Western blotting or mass spectrometry (MS).



**Figure 27. Easy RNA labeling and interaction analysis.** RNA probes are first labeled at the 3' end with desthiobiotin using T4 RNA Ligase. RNA probes are then immobilized onto magnetic streptavidin particles and incubated with protein from cell lysates or *in vitro* translation preps. RNA-binding proteins are eluted and detected by Western blotting or mass spec.

### Highlights

- **Direct** – capture ribonucleoprotein complexes directly from cell lysates
- **Easy to use** – magnetic beads enable easy processing for multiple samples
- **Flexible** – enrich RNA binding proteins from cell/tissue lysates or *in vitro* translated protein preps
- **Clean** – magnetic format yields low background
- **Specific** – perform RNA mutations to map interaction sites
- **Complete** – contains both labeling and enrichment modules with buffers necessary for assay; positive control RNA, negative control RNA and anti-RBP antibody included



**Figure 28. End-labeled RNA enriches specific target binding proteins.** RNA binding proteins (RBP) of the AR 3' UTR control system (top panel) and three experimental systems were enriched according to kit procedure. L = lysate; FT = flow-through; E = elute.

Incubation of A431 lysate with labeled AR UTR RNA (Target) enriches HuR RBP, while incubation with negative control poly(A) RNA (Unrelated) or beads only (None) does not (compare elution lanes). The same pattern results with the experimental systems, confirming the proper function of the kit. Samples were normalized by volume, and bands were detected using Thermo Scientific™ SuperSignal™ West Pico Substrate by a 2-minute exposure to film. Target RNA sequences were as follows:

Androgen Receptor 3' UTR (Kit Control System):  
 5'-CUGGGCCUUUUUUUCUCUCCUCCUUCUUUUUUUCUUCUCCUCCUA-3'  
 Tat Exon 2:  
 5'-UUACUCAACAGAGGAGAGCAAGAAUUGGAGCCAGUAGAUCCUAGACUAGAGCCUUGG-3'  
 Let-7 loop:  
 5'-CAGUUUGAGGGUCUAUGAUACCACCGGUACAAGAUACUG-3'  
 Poly(A)25:  
 5'-AAAAAAAAAAAAAAAAAAAAAAAAA-3'



# perform **high-capacity** protein purification


## Magnetic Biotin Pull-Down

### Thermo Scientific™ Pierce™ Streptavidin Magnetic Beads



provide easy affinity purification of biotin-labeled target molecules without columns or centrifugation. Pierce Streptavidin Magnetic Beads use a recombinant form of streptavidin with a mass of 53kDa and a near-neutral isoelectric point (pI). The protein is a tetramer having four biotin-binding sites. Unlike avidin, streptavidin has no carbohydrate groups, resulting in low nonspecific binding. The high-affinity interaction between streptavidin and biotin cannot be dissociated efficiently except under very harsh conditions, such as boiling in sample loading buffer for SDS-PAGE or 8 M guanidine•HCl, pH 1.5. Consequently, it is often possible to elute binding partners in an interaction complex without also eluting the biotinylated component.

#### Highlights

- **Stable immobilization chemistry** – streptavidin is immobilized using leach-resistant chemistry
- **High capacity** – superior-quality beads with high binding capacity provide rapid and efficient biomolecule purification from complex samples
- **Low nonspecific binding** – stable, pre-blocked beads provide clean purification products that are compatible with mass spectrometry analysis
- **Superior performance** – binding capacity is nearly three times higher than beads from other suppliers, allowing the use of smaller samples per experiment



**Figure 31. Better immunoprecipitation results with Thermo Scientific Pierce Streptavidin Magnetic Beads.** MOPC cell lysate (0.75mg per sample) was incubated overnight at 4°C with and without 10µg biotinylated Grp94 antibody. Pierce Streptavidin Magnetic Beads and Dynabeads MyOne™ Streptavidin T1 Beads (Life Technologies) were added to a 96 deep-well plate (0.5mg or 0.25mg per well). Eluates were resolved by SDS-PAGE and analyzed by Western blot with anti-Grp94 antibody. About 0.25mg of Pierce Streptavidin Magnetic Beads gave the same yield as 0.5mg of Dynabeads MyOne T1 Beads.



**Figure 32. Higher binding capacity with Thermo Scientific Pierce Streptavidin Magnetic Beads.** Pierce Streptavidin Magnetic Beads and Life Technologies Dynabeads MyOne Streptavidin T1 Beads were added to a 96 deep-well plate (1mg beads per well). Using the KingFisher 96 Instrument, the beads were washed with phosphate-buffered saline containing 0.05% Tween-20. The beads were then incubated for one hour with varying amounts of biotinylated rabbit IgG (20-225µg).



#### Ordering Information



Product #	Description	Pkg. Size
88816	<b>Pierce Streptavidin Magnetic Beads</b> Sufficient for: Binding approx. 55µg biotinylated rabbit IgG per mg of beads (approx. 3500 pmol biotinylated fluorescein per mg of beads)	1 mL
88817	<b>Pierce Streptavidin Magnetic Beads</b> Sufficient for: Binding approx. 55µg biotinylated rabbit IgG per mg of beads (approx. 3500 pmol biotinylated fluorescein per mg of beads)	5 mL

## Magnetic GST-Tagged Protein Purification

**Thermo Scientific™ Pierce™ Glutathione Magnetic Beads** provide a simple, rapid and reliable method for the purification of glutathione S-transferase (GST) fusion proteins from crude cell lysate prepared from bacteria, yeast or mammalian cells. These beads can be used to isolate GST-tagged proteins or perform pull-down assays using GST-tagged proteins as bait.

### Highlights

- **High binding** – 5-10mg GST/mL settled beads
- **Stable affinity ligand** – glutathione is covalently immobilized to particles, ensuring clean, leach-resistant purification products
- **High capacity** – binding capacity is sufficient for both routine and demanding magnetic separation procedures



**Figure 33. High-performance purification of a GST fusion protein using Thermo Scientific Pierce Glutathione Magnetic Beads.** Bacterial cells expressing GST-Rabaptin were lysed, and replicate aliquots were processed with the Pierce Glutathione Magnetic Beads in a 96 deep-well plate using a KingFisher 96 Instrument. Eluates were boiled in reducing sample buffer, resolved by SDS-PAGE and stained with coomassie dye. Purity and reproducibility were excellent.



### Ordering Information

Product #	Description	Pkg. Size
88821	<b>Pierce Glutathione Magnetic Beads</b> <i>Sufficient for: Binding 5 to 10mg GST per mL of beads</i>	4mL
88822	<b>Pierce Glutathione Magnetic Beads</b> <i>Sufficient for: Binding 5 to 10mg GST per mL of beads</i>	20mL


# perform **high-capacity** protein purification

## Magnetic His-Tagged Protein Purification

**Thermo Scientific™ HisPur™ Ni-NTA Magnetic Beads** are high-capacity Nickel-IMAC beads for affinity purification of His-tagged fusion proteins in manual or automated formats. The blocked magnetic bead surface is derivatized with the nitrilotriacetic acid (NTA) chelation moiety and loaded with divalent nickel ions ( $\text{Ni}^{2+}$ ). The immobilized metal affinity chromatography (IMAC) beads provide high binding capacity with very low background. The HisPur Ni-NTA Magnetic Beads can be used both manually with a magnetic stand as well as with automated platforms such as the KingFisher Instruments for high-throughput needs.

### Highlights



- **High capacity** – equivalent or higher binding capacity than Ni-NTA magnetic beads from other suppliers
- **Low nonspecific binding** – the bead surface is pre-blocked and the protocol provides optimized buffers for purification
- **Fast** – protocol is completed in less than one hour
- **Scalable** – process microliter to milliliter sample volumes
- **Versatile** – purify proteins using native or denaturing conditions
- **Reagent compatible** – can be used with common cell lysis reagents and a variety of buffer additives
- **Multiple formats** – protein coupling to the beads and downstream applications can be performed both manually and on an automated platform (e.g., KingFisher Instruments)



**Figure 34. Thermo Scientific HisPur Ni-NTA Magnetic resin delivers consistent yield.**



His-tag protein purification was performed in a 96-well plate using a KingFisher Flex Instrument. In each well, 100µg of *E. coli* lysate expressing 6xHis-GFP protein was added to 0.5mg of Thermo Scientific™ Pierce™ Magnetic Ni-NTA Resin. Eluted protein was analyzed by SDS-PAGE stained with Imperial Protein Stain to determine well-to-well consistency in protein recovery. The variance between samples is measured at less than 15%.





**Table 8. Characteristics of Thermo Scientific HisPur Ni-NTA Magnetic Beads.**

<b>Composition</b>	Ni <sup>2+</sup> loaded on nitrilotriacetic acid that has been covalently coupled to the beads
<b>Mean Diameter</b>	1µm (nominal)
<b>Density</b>	2.0g/cm <sup>3</sup>
<b>Bead Concentration</b>	12.5mg/mL in 20% ethanol
<b>Binding Capacity</b>	≥40µg of 6X-His-tagged GFP/mg of beads; ≥500µg of 6X-His-tagged GFP/mL of beads



### Ordering Information

Product #	Description	Pkg. Size
88831	HisPur Ni-NTA Magnetic Beads	2mL
88832	HisPur Ni-NTA Magnetic Beads	10mL




**Figure 35. Superior performance of Thermo Scientific HisPur Ni-NTA Magnetic Beads.** Bacterial lysate (100µL total protein) containing over-expressed 6xHis-GFP (**Panel A**) or over-expressed 6xHis-β Galactosidase (**Panel B**) was applied to 0.5mg of HisPur Ni-NTA Magnetic Beads, Mag Sepharose (GE Life Sciences), PureProteome (EMD/Millipore) or Magnetic Agarose (Qiagen) Ni-NTA beads. All samples were run in duplicate, and the beads were processed using buffers recommended by the manufacturers. For the HisPur Ni-NTA Magnetic Beads, the amount of imidazole in the equilibration, wash and elution buffers was 30mM, 50mM and 150mM, respectively. All three buffers contained 100mM sodium phosphate and 600mM sodium chloride. Binding was performed with all samples for 30 minutes. The beads were collected on a magnetic stand and the flow-throughs were saved for analysis. Eluates were resolved on an SDS-PAGE gel and stained with Imperial Stain. For purification of 6xHis-GFP, comparable yields and purity were observed for HisPur Ni-NTA and Qiagen Ni-NTA Magnetic Beads. HisPur Ni-NTA Magnetic Beads showed higher yield and purity than the Qiagen Magnetic Agarose Ni-NTA beads in the purification of 6xHis-β Galactosidase. The Mag Sepharose and PureProteome Ni-NTA Magnetic beads gave lower purity and lower yield than HisPur Ni-NTA Magnetic Beads in both purifications. Purity analyses were performed on a Thermo Scientific™ mECL™ Imager with Thermo Scientific™ mImageAnalysis™ Software. Purity was determined by measuring the ratio of the background-corrected 6xHis-tagged protein band of interest to the sum of all bands in a given lane.

# isolate phosphopeptides for mass spec analysis

## Magnetic Phosphopeptide Enrichment

### The Thermo Scientific™ Pierce™ Magnetic Titanium Dioxide Phosphopeptide Enrichment Kit

isolates phosphopeptides from complex biological samples using titanium dioxide-coated magnetic beads. The TiO<sub>2</sub> ligand selectively binds peptides containing phosphorylated serine (Ser), tyrosine (Tyr) or threonine (Thr), enabling phosphopeptide enrichment from protease-digested samples.



**Table 9. Phosphopeptide enrichment improves MS identification of phosphoproteins.** Two milligrams of a tryptic digest prepared from peripheral blood mononuclear cells (lymphocytes) with and without phosphopeptide enrichment were analyzed by MS. Enrichment was performed with the Pierce Titanium Dioxide Phosphopeptide Enrichment Kit using the KingFisher 96 Instrument. Samples were analyzed on a Thermo Scientific™ LTQ Orbitrap™ Mass Spectrometer.

	Enriched	Non-Enriched
Total number of proteins identified	185	247
Total number of phosphoproteins identified	160	1
Total number of peptides identified	2347	2457
Total number of phosphopeptides identified	2009	7
Total number of unique phosphopeptides identified	177	1
Relative enrichment for phosphopeptides (%)	86	0.3

The isolated phosphopeptides are compatible for analysis downstream by mass spectrometry (MS).

#### Highlights

- **Complete MS-compatible kits** – include ready-to-use binding, wash and elution buffers that are optimized for phosphopeptide enrichment and downstream analysis by MALDI- and ESI-based MS
- **Optimized for high-throughput screening (HTS)** – validated procedure for processing from 1 to 96 samples at a time; complete entire assay in about 15 minutes using a KingFisher Flex Instrument
- **Stable affinity ligand** – titanium dioxide is specially coated as a film onto the magnetic particles
- **Selective** – affinity system is selective for phosphorylated Ser, Tyr and Thr; exhibits minimal nonspecific binding to acidic residues
- **Sensitive** – affinity provides more than 1000 times greater sensitivity than traditional IMAC technologies; enables enrichment and MS measurement of less than 100fmol of phosphoprotein



#### Ordering Information

Product #	Description	Pkg. Size
88811	<b>Pierce Magnetic Titanium Dioxide Phosphopeptide Enrichment Kit</b> <i>Sufficient for: Purifying 96 x 100µg peptide samples</i> Contains: TiO <sub>2</sub> Magnetic Beads (20X), 1mL Binding Buffer, 100mL Washing Buffer, 25mL Elution Buffer, 3mL Thermo-Fast 96 Robotic PCR Plates, 0.2mL wells, 2 plates	96-rxn kit
88812	<b>Pierce Magnetic Titanium Dioxide Phosphopeptide Enrichment Kit, Trial Size</b> <i>Sufficient for: Purifying 24 x 100µg peptide samples</i> Contains: TiO <sub>2</sub> Magnetic Beads (20X), 0.25mL Binding Buffer, 100mL Washing Buffer, 25mL Elution Buffer, 3mL Thermo-Fast 96 Robotic PCR Plates, 0.2mL wells, 2 plates	24-rxn kit

# high-throughput quantitative proteomics

## Magnetic Bead Processors and Mass Spectrometers

**Our revolutionary proprietary magnetic separation technology** lets you process virtually any sample from any source for the ultimate in isolation of nucleic acids, proteins and cells. With four platforms to choose from, the Thermo Scientific™ KingFisher Systems provide the performance, flexibility and speed for your budget, application and throughput requirements.



### Thermo Scientific KingFisher Flex System

With high-throughput or processing volume of up to 5mL, the KingFisher Flex System offers truly versatile purification of nucleic acids and proteins. Process volumes from 20µL to 5000µL,

depending on the magnet head, with 96- or 24-well format. Use predefined protocols or customize your own for special applications. The new KingFisher Flex System replaces the KingFisher 96 Instrument.

### Thermo Scientific KingFisher Duo System

New to the KingFisher family, the KingFisher Duo System delivers advanced functionality in a compact, mid-throughput capacity instrument for isolation applications. Its small footprint and big functionality, including traceability and data management, make it a perfect fit for research and routine laboratories. Two protocols can run sequentially without interruption, raising throughput up to 24 samples per load. The KingFisher Duo System also includes large volume processing of up to 5mL.



### Thermo Scientific KingFisher mL System

The Thermo Scientific KingFisher mL System enables automated, low-throughput sample preparation into your laboratory workflow. Processes volumes from 50µL to 1000µL carried out using tube strips.



### Original Thermo Scientific KingFisher System

The first in the family, the Thermo Scientific™ KingFisher™ System allows you to economically purify small-scale samples. Run up to 24 samples of 20µL to 200µL. All purification and processing steps can be programmed using simple push-button operation and are carried out in microstrips.

### Thermo Scientific KingFisher Kits

With optimized Thermo Scientific™ KingFisher™ Purification Kits, you can easily perform blood DNA, total RNA, cell and tissue DNA, viral NA, and plant DNA extraction. KingFisher Instruments, Software, Kits and Consumables deliver unparalleled performance.



### Thermo Scientific Orbitrap and Orbitrap Hybrid Mass Spectrometers

- Thermo Scientific™ Exactive™ Plus MS
- Thermo Scientific™ Q Exactive™ Hybrid Quadrupole-Orbitrap MS
- Thermo Scientific™ LTQ Orbitrap XL™ Hybrid Ion Trap-Orbitrap MS
- Thermo Scientific™ Orbitrap™ Velos Pro Hybrid Ion Trap-Orbitrap MS
- Thermo Scientific™ Orbitrap Elite Hybrid Ion Trap-Orbitrap MS

Instrument	Flex		Duo		mL	KF
Samples/run	96	24	12(24)	6	15	24
Working volume (µL)	20-1000	200-5000	30-1000	200-5000	50-1000	20-200

For more information on Thermo Scientific KingFisher Systems, visit [thermoscientific.com/kingfisher](http://thermoscientific.com/kingfisher) or consult your local sales representative.

For more information on Thermo Scientific Mass Spectrometers, visit [thermoscientific.com/ms](http://thermoscientific.com/ms) or consult your local sales representative.


# your helping hands for protein research

To view a PDF or request a hard copy of one of these handbooks, visit [thermoscientific.com/pierce](http://thermoscientific.com/pierce) and click on Technical Resources > Request Literature.




## Cell Lysis Technical Handbook

The Cell Lysis Technical Handbook describes the latest Thermo Scientific Cell Lysis products. Included are novel lysis products for neuronal cells and synaptosomes, subcellular protein fractionation kits from tissues, new protease and/or phosphatase inhibitor tablets, and universal nuclease to reduce sample viscosity.



## Mass Spec Sample Preparation Handbook

This handbook provides background, helpful hints and troubleshooting advice for cell lysis, sample preparation, detection, mass spectrometry sample preparation and downstream applications. The handbook features new products for protein concentration, purification and enrichment, plus the latest labeling techniques, including Thermo Scientific™ SILAC, TMT™, cystTMT™ and HeavyPeptide™ Reagents. The book also includes a section on Thermo Scientific Mass Spectrometry Instrumentation and Software. Everything you need to extract, digest, enrich, clean up and quantify proteins and peptides in one volume.



## Western Blotting Handbook and Troubleshooting Guide

The Western Blotting Handbook and Troubleshooting Guide (version 3) details each step of the Western blotting process with technical information and products for transfer, blocking, washing, antibodies, substrates, film and stripping buffer. You will want to keep this booklet close at hand because it also includes protocols, references and a troubleshooting guide.



Find us on  
Facebook

[thermoscientific.com/pierce](http://thermoscientific.com/pierce)

© 2013 Thermo Fisher Scientific Inc. All rights reserved. These products are supplied for laboratory or manufacturing applications only. Dynabeads and MyOne are trademarks or registered trademarks of Life Technologies Corp. Sepharose is a trademark of GE Healthcare Bio-Sciences. PureProteome is a trademark of Merck KGaA. Tween is a trademark of ICI Americas. Tandem Mass Tag, TMT and cystTMT are trademarks of Proteome Sciences plc. Mascot is a trademark of Matrix Science Ltd. SPHERO is a trademark of Spherotech, Inc. IQ5 is a trademark of Bio-Rad Laboratories, Inc. ProteoPep is a trademark of New Objective, Inc. Facebook is a trademark of Facebook, Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

### Life Science Research

**Africa/Belgium/Europe/Middle East**

+32 53 85 71 84

**France** 0 800 50 82 15

**Germany** 0228 9125650

**Netherlands** 076 50 31 880

**Switzerland** 0800 56 31 40

**UK** 0800 252 185

Email: [perbio.euomarketing@thermofisher.com](mailto:perbio.euomarketing@thermofisher.com)  
[thermoscientific.com/pierce](http://thermoscientific.com/pierce)

For other regions, visit  
[thermoscientific.com/piercedistributors](http://thermoscientific.com/piercedistributors)

**USA** +815-968-0747 or +800-874-3723

Customer Assistance

E-mail: [Pierce.CS@thermofisher.com](mailto:Pierce.CS@thermofisher.com)  
[thermoscientific.com/pierce](http://thermoscientific.com/pierce)

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific