

GE Healthcare

# His SpinTrap

50 prepacked His SpinTrap columns

Product booklet

Codes: 28-4013-53  
28-9321-71



# Page finder

|   |    |
|---|----|
| 1. Legal  | 3  |
| 2. Handling   | 4  |
| 2.1. Safety warnings and precautions                                | 4  |
| 2.2. Storage  | 4  |
| 3. Purpose  | 5  |
| 4. Principle  | 6  |
| 5. Advice on handling   | 9  |
| 5.1. Buffers for native conditions                                  | 9  |
| 5.2. Buffers for denaturing conditions                              | 9  |
| 5.3. Optimal imidazole concentration in binding and elution buffers | 9  |
| 6. Sample pretreatment protocol                                     | 10 |
| 6.1. Dilute the cell paste  | 10 |
| 6.2. Enzymatic lysis  | 10 |
| 6.3. Mechanical lysis   | 10 |
| 6.4. Clarify the lysate   | 11 |
| 7. Purification protocol  | 12 |
| 8. Prepare buffers with different imidazole concentrations          | 14 |
| 8.1. 2 M imidazole stock solution                                   | 14 |
| 8.2. Eluted histidine-tagged protein not pure                       | 14 |
| 9. Tips and hints   | 16 |
| 9.1. Liquid not completely removed during centrifugation            | 16 |
| 9.2. Eluted histidine-tagged protein not pure                       | 16 |
| 9.3. Low yield of histidine-tagged protein                          | 17 |
| 10 Ordering information   | 18 |

# 1. Legal

GE, imagination at work and GE monogram are trademarks of General Electric Company.

SpinTrap, GraviTrap, MultiTrap, HisTrap, and Sepharose are trademarks of GE Healthcare companies.

All third party trademarks are the property of their respective owners.

This product is covered by US pat No 6 623 655 and their equivalents in other countries.

© 2005-2010 General Electric Company – All rights reserved.  
Previously published Sep. 2005

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request.

Contact your local GE Healthcare representative for the most current information.

<http://www.gelifesciences.com/sampleprep>

GE Healthcare UK Limited.  
Amersham Place, Little Chalfont,  
Buckinghamshire, HP7 9NA UK

## 2. Handling

### 2.1. Safety warnings and precautions

**Warning: For research use only.**

Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

All chemicals should be considered as potentially hazardous. Only persons trained in laboratory techniques and familiar with the principles of good laboratory practice should handle this product. Suitable protective clothing such as laboratory overalls, safety glasses, and gloves should be worn. Care should be taken to avoid contact with skin or eyes; if contact should occur, wash immediately with water. See material safety data sheet(s) and/or safety statement(s) for specific advice.

**Warning: .**

The column contains nickel, Ni<sup>2+</sup>, which is potentially allergenic. Always use normal personal protection devices like gloves and safety glasses when handling His SpinTrap columns.

**Warning: .**

The column storage solution, 0.15% Kathon CG/ICP Biocide, is potentially allergenic. Use gloves when discarding the storage solution.

### 2.2. Storage

Store at 4°C to 30°C.

### 3. Purpose

His SpinTrap™ columns are designed for rapid small-scale purification of histidine-tagged proteins. The columns are suitable for purification of multiple samples in parallel, for example in screening experiments. For increased convenience a combination kit is also available consisting of His SpinTrap and His Buffer kit.

**His SpinTrap contains:**

- 50 prepacked His SpinTrap™ columns
- Instructions for use

## 4. Principle

His SpinTrap contains Ni Sepharose™ High Performance medium, which has high binding capacity for histidine-tagged proteins. Table 1 summarizes His SpinTrap characteristics.

**Table 1.** His SpinTrap characteristics.

|                                       |   |
|---------------------------------------|---|
| Column material                       | Polypropylene barrel, polyethylene frits  |
| Medium                                | Ni Sepharose High Performance   |
| Average bead size                     | 34 µm   |
| Protein binding capacity <sup>1</sup> | Approx. 750 µg histidine-tagged protein/column  |
| Bed volume                            | 100 µl  |
| Compatibility during use              | Stable in all commonly used buffers, reducing agents, denaturants, and detergents (See Table 2) |
| Avoid in buffers                      | Chelating agents, e.g. EDTA, EGTA, citrate. (See Table 2)                                       |
| Storage                               | 0.15% Kathon™ CG  |
| Storage temperature                   | 4°C to 30°C   |

<sup>1</sup> Binding capacity is protein-dependent.

The column has low nickel ion (Ni<sup>2+</sup>) leakage and is compatible with denaturing agents and a wide range of additives, see Table 2.

**Table 2.** His SpinTrap compatibility.

|                                |   |
|--------------------------------|---|
| Reducing agents                | 5 mM DTE                                  |
|                                | 5 mM DTT                                  |
|                                | 20 mM $\beta$ -mercaptoethanol            |
|                                | 5 mM TCEP                                 |
|                                | 10 mM reduced glutathione                 |
| Denaturing agents <sup>1</sup> | 8 M urea                                  |
|                                | 6 M guanidine-HCl                         |
| Detergents                     | 2% Triton™ X-100 (nonionic)               |
|                                | 2% Tween™ 20 (nonionic)                   |
|                                | 2% NP-40 (nonionic)                       |
|                                | 2% cholate (anionic)                      |
|                                | 1% CHAPS (zwitterionic)                   |
| Other additives                | 20% ethanol                               |
|                                | 50% glycerol                              |
|                                | 100 mM Na <sub>2</sub> SO <sub>4</sub>    |
|                                | 1.5 M NaCl                                |
|                                | 1 mM EDTA <sup>2</sup>                    |
|                                | 60 mM citrate <sup>2</sup>                |
| Buffers                        | 50 mM sodium phosphate, pH 7.4            |
|                                | 100 mM Tris-HCl, pH 7.4                   |
|                                | 100 mM Tris-acetate, pH 7.4               |
|                                | 100 mM HEPES, pH 7.4                      |
|                                | 100 mM MOPS, pH 7.4                       |
|                                | 100 mM sodium acetate <sup>1</sup> , pH 4 |

<sup>1</sup> Tested for one week at 40°C.

<sup>2</sup> Generally, chelating agents should be used with caution (and only in the sample, not in the buffers). Any metal-ion stripping may be counteracted by adding a small excess of MgCl<sub>2</sub> before centrifugation/filtration of the sample.

His SpinTrap can be used with a standard microcentrifuge, and a purification takes approximately 10 minutes.

Cell culture lysates may be directly applied to the column without prior clarification.



## 5. Advice on handling

### 5.1. Buffers for native conditions

Recommended buffers for native conditions can easily be prepared from His Buffer Kit or according to the description in Section 6.

- Binding buffer  
20 mM sodium phosphate, 500 mM NaCl, 20 mM imidazole, pH 7.4
- Elution buffer  
20 mM sodium phosphate, 500 mM NaCl, 500 mM imidazole, pH 7.4

### 5.2. Buffers for denaturing conditions

- Binding buffer  
20 mM Tris-HCl, 8 M urea, 500 mM NaCl, 5 mM imidazole, pH 8.0  
+ 1 mM  $\beta$ -mercaptoethanol
- Elution buffer  
20 mM Tris-HCl, 8 M urea, 500 mM NaCl, 500 mM imidazole, pH 8.0  
+ 1 mM  $\beta$ -mercaptoethanol

### 5.3. Optimal imidazole concentration in binding and elution buffers

#### **Binding buffer**

The optimal imidazole concentration in the binding buffer is protein dependent and has influence on final yield and purity of the histidine-tagged protein. Under native conditions, 20–40 mM imidazole in the binding buffer is suitable for many proteins.

#### **Elution buffer**

A concentration of 500 mM imidazole in the elution buffer is most often sufficient to completely elute the target protein.

**Note:** As an alternative to elution with imidazole, pH can be lowered to approximately pH 4.5. At pH below 4.0 the metal ions will be stripped off the medium.

## 6. Sample pretreatment protocol

This is the recommended sample pretreatment protocol. However, other established sample pretreatment procedures may also work.

Use standard 2 ml microcentrifuge tubes.

### 6.1. Dilute the cell paste

- Add 1 ml binding buffer to resuspend cell paste obtained from 20 to 50 ml cell culture (volume depending on expression level).

**Note:** To prevent binding of host cell proteins it is essential that the sample and binding buffers contain the same concentration of imidazole.

### 6.2. Enzymatic lysis

- Add the following substances to specified final concentrations in the cell suspensions:
  - Lysozyme: 0.2 mg/ml
  - DNase: 20 µg/ml
  - MgCl<sub>2</sub>: 1 mM
  - Pefabloc™ SC or PMSF: 1 mM
- Vortex the tubes gently and incubate at room temperature for 30 min.

**Note:** Chemical lysis kits can also be used, but make sure that they do not contain any chelating agent.

### 6.3. Mechanical lysis

- Repeated freeze/thaw or sonication.

## 6.4. Clarify the lysate

- Spin at full speed in a microcentrifuge for 10 min to remove insoluble material.
- Collect supernatants and purify on His SpinTrap.

**Note:** Cell culture lysates may be directly applied to the column without prior clarification (i.e. omit step 4).

## 7. Purification protocol

Run purifications on His SpinTrap using a standard microcentrifuge. Place the column in a 2 ml microcentrifuge tube to collect the liquid during centrifugation.

Use a new 2 ml tube for every step.

---

### 1. Remove storage solution

- Invert and shake the column repeatedly to resuspend the medium.
- Loosen the top cap one-quarter of a turn and twist off the bottom closure.
- Place the column in a 2 ml microcentrifuge tube and centrifuge for 30 s at 70 to 100  $\times$  g.
- Remove and discard the top cap.



---

### 2. Column equilibration

- Add 600  $\mu$ l binding buffer.
- Centrifuge for 30 s at 70 to 100  $\times$  g




---

### 3. Sample application

- Add up to 600  $\mu$ l sample in one application.
- Centrifuge for 30 s at 70 to 100  $\times$  g.

**Note:** Several sample applications can be performed as long as the capacity of the column is not exceeded.

---



---

#### 4. Wash

- Add 600  $\mu$ l binding buffer.
- Centrifuge for 30 s at 70 to 100  $\times$  g.



---

#### 5. Elution

- Add 200  $\mu$ l elution buffer.
- Centrifuge for 30 s at 70 to 100  $\times$  g and collect the purified sample.
- Add 200  $\mu$ l elution buffer.
- Centrifuge for 30 s at 70 to 100  $\times$  g and collect the purified sample.



**Note:** The first eluted 200  $\mu$ l will contain the majority of the target protein.

---

## 8. Prepare buffers with different imidazole concentrations

### 8.1. 2 M imidazole stock solution

To prepare binding and elution buffer with different amounts of imidazole a 2 M imidazole stock solution should be used.

Prepare 250 ml 2 M imidazole stock solution:

1. To 34.05 g imidazole, add distilled water to 200 ml and dissolve completely. Use high purity imidazole as this will give no or very low absorbance at 280 nm (imidazole, 68.08 g/mol).
2. Adjust to pH 7.4 with HCl.
3. Add distilled water to 250 ml.

### 8.2. Binding or elution buffers

Prepare 250 ml binding or elution buffers with final concentrations, 20 mM sodium phosphate, 500 mM NaCl and 10–500 mM imidazole:

1. Take 0.44 g  $\text{Na}_2\text{HPO}_4 \times 2\text{H}_2\text{O}$  (177.99 g/mol), 0.35 g  $\text{NaH}_2\text{PO}_4 \times \text{H}_2\text{O}$  (137.99 g/mol) and 7.30 g NaCl (58.44 g/mol).
2. Add X ml 2 M imidazole stock solution (see Table 3). The volume of imidazole stock solution added depends on the chosen imidazole binding and elution concentrations.
3. Add distilled water to 200 ml and dissolve completely.
4. Adjust to pH 7.4 with HCl.
5. Add distilled water to 250 ml.
6. Filter through a 0.45  $\mu\text{m}$  filter.

**Table 3.** Phosphate buffers (250 ml) with different imidazole concentrations.

| <b>Final imidazole concentration</b> | <b>2M imidazole stock solution volume (ml)</b> |
|--------------------------------------|--|
| 10                                   | 1.25   |
| 20                                   | 2.50   |
| 30                                   | 3.75   |
| 40                                   | 5.00   |
| 50                                   | 6.25   |
| 60                                   | 7.50   |
| 70                                   | 8.75   |
| 80                                   | 10.00  |
| 90                                   | 11.25  |
| 100                                  | 12.50  |
| 300                                  | 37.50  |
| 500                                  | 62.50  |

## 9. Tips and hints

### 9.1. Liquid not completely removed during centrifugation

#### **Sample too viscous**

- Increase centrifugation time.
- Increase dilution of the cell paste before or after mechanical lysis.
- Continue mechanical lysis until the viscosity is reduced, and/or add an additional dose of DNase and Mg<sup>2+</sup>.
- Filter the sample (or centrifuge if you have used unclarified sample).

#### **Target protein difficult to dissolve or precipitates**

- Add detergents, reducing agents or other additives (Table 2) and mix gently for 30 min to aid solubilization of the tagged protein. Note that Triton X-100 and NP-40 (but not Tween) have a high absorbance at 280 nm. Furthermore, detergents cannot be easily removed by buffer exchange.
- Inclusion bodies: the protein can usually be solubilized (and unfolded) from inclusion bodies using common denaturants such as 4–6 M guanidine-HCl, 4–8 M urea or strong detergents. Mix gently for 30 min or more to aid solubilization.

### 9.2. Eluted histidine-tagged protein not pure

#### **Imidazole concentration in sample and binding buffer too low**

- Increase imidazole concentration in sample and binding buffer to prevent contaminants binding. We recommend 20–40 mM, but higher concentrations may also work.

#### **Partial degradation of tagged protein by proteases**

- Add protease inhibitors (use EDTA with caution, see Table 2). Perform lysis and purify at 4°C.



### **Contaminants are associated with tagged proteins**

- Add detergent and/or reducing agents before sonicating the cells. Increase detergent levels (e.g. up to 2% Triton X-100 or 2% Tween), or add glycerol (up to 50%) to the wash buffer to disrupt nonspecific interactions.

### **Insufficient washing of unbound material**

- Repeat the wash step after sample application to obtain optimal purity.

## **9.3. Low yield of histidine-tagged protein**

### **Histidine-tagged protein found in the flowthrough during sample application and wash**

- Imidazole concentration in the sample and binding buffer is too high. Use a lower concentration.
- Ensure that the concentration of chelating or strong reducing agents in the sample is not too high.
- The histidine tag may be insufficiently exposed; perform purification of unfolded protein in urea or guanidine-HCl as for inclusion bodies. To minimize dilution of the sample, add solid urea or guanidine-HCl.
- The histidine tag has been lost. Check the sequence of the construct.

### **Histidine-tagged protein not eluted during purification**

- Histidine-tagged protein still bound. Elute with a higher concentration of imidazole in the elution buffer.
- The target protein has precipitated in the column. Decrease the amount of sample. Decrease imidazole concentration during elution. Try detergents or change NaCl concentration, or elute under denaturing (unfolding) conditions.
- Nonspecific hydrophobic or other interaction. Add a nonionic detergent to the elution buffer or increase NaCl concentration.

### **Histidine-tagged protein not completely eluted**

- Elute with a larger volume of elution buffer.

## 10. Ordering information

| <b>Product</b>                | <b>Pack size</b> | <b>Code No.</b> |
|-------------------------------|------------------|-----------------|
| His SpinTrap                  | 50 × 100 µl      | 28-4013-53      |
| His SpinTrap Kit <sup>1</sup> | 1                | 28-9321-71      |

| <b>Related products</b>          | <b>Pack size</b>                | <b>Code No.</b> |
|----------------------------------|---------------------------------|-----------------|
| His Buffer Kit                   | 1                               | 11-0034-00      |
| His GraviTrap™                   | 10 × 1 ml                       | 11-0033-99      |
| His GraviTrap Kit <sup>2</sup>   | 1                               | 28-4013-51      |
| His MultiTrap™ FF                | 4 × prepacked<br>96-well plates | 28-4009-90      |
| His MultiTrap HP                 | 4 × prepacked<br>96-well plates | 28-4009-89      |
| HisTrap™ FF, 1 ml                | 5 × 1 ml                        | 17-5319-01      |
| HisTrap FF, 5 ml                 | 5 × 5 ml                        | 17-5255-01      |
| HisTrap FF crude, 1 ml           | 5 × 1 ml                        | 11-0004-58      |
| HisTrap FF crude, 5 ml           | 5 × 5 ml                        | 17-5286-01      |
| HisTrap HP, 1 ml                 | 5 × 1 ml                        | 17-5247-01      |
| HisTrap HP, 5 ml                 | 1 × 5 ml                        | 17-5248-01      |
| HisTrap HP, 5 ml                 | 5 × 5 ml                        | 17-5248-02      |
| Ni Sepharose 6 Fast Flow         | 5 ml                            | 17-5318-06      |
| Ni Sepharose 6 Fast Flow         | 25 ml                           | 17-5318-01      |
| Ni Sepharose 6 Fast Flow         | 100 ml                          | 17-5318-02      |
| Ni Sepharose High<br>Performance | 25 ml                           | 17-5268-01      |

| <b>Related products</b>        | <b>Pack size</b> | <b>Code No.</b> |
|--------------------------------|------------------|-----------------|
| Ni Sepharose High Performance  | 100 ml           | 17-5268-02      |
| Empty Disposable PD-10 columns | 50               | 17-0435-01      |
| LabMate PD-10 Buffer Reservoir | 10               | 18-3216-03      |

<sup>1</sup> Includes 1 pack His SpinTrap and 1 pack His Buffer Kit

<sup>2</sup> Includes 2 packs His GraviTrap and 1 pack His Buffer Kit

| <b>Literature</b>  | <b>Code No.</b> |
|--|-----------------|
| Ni Sepharose and IMAC Sepharose, Selection guide           | 28-4070-92      |
| Affinity Chromatography Columns and Media, Selection guide | 18-1121-86      |
| Affinity Chromatography Handbook, Principle and Methods    | 18-1022-29      |
| Recombinant Protein Purification Handbook                  | 18-1142-75      |
| Data File His SpinTrap and His SpinTrap Kit                | 28-4046-59      |

**GE Healthcare offices:**

GE Healthcare Bio-Sciences AB  
Björkgatan 30, 751 84 Uppsala,  
Sweden

GE Healthcare Europe GmbH  
Munzinger Strasse 5, D-79111 Freiburg,  
Germany

GE Healthcare Bio-Sciences Corp.  
800 Centennial Avenue, P.O. Box 1327,  
Piscataway, NJ 08855-1327,  
USA

GE Healthcare Japan Corporation  
Sanken Bldg. 3-25-1, Hyakunincho,  
Shinjuku-ku, Tokyo 169-0073,  
Japan

For contact information for your local office,  
please visit: [www.gelifesciences.com/contact](http://www.gelifesciences.com/contact)

GE Healthcare UK Limited  
Amersham Place  
Little Chalfont, Buckinghamshire,  
HP7 9NA, UK

<http://www.gelifesciences.com/sampleprep>



imagination at work