thermoscientific



Thermo Scientific Samco Clicktainer Vial

Unique click seal system providing sample protection and user safety





Thermo Scientific™ Samco™ Clicktainer™ Vial

Providing ultimate sample protection and user safety by ensuring correct cap application the first time

Thermo Fisher Scientific understands the value of patient samples. Sample collection needs to be easy. Sample transport needs to be secure. The Samco Clicktainer vials have a wide opening and leakproof seal. The "click" seal system indicates when the cap is properly applied so the sample is secure.





User Safety and Sample Protection



Clicktainer cap

- Visual, audible and tactile indication of full cap closure
- Safe tab design reduces risk of tear to gloves or break off during use

Pneumatic Transport Tested

 Ideal for use in hospital pneumatic transport systems

Excellent leakproof performance

- Pressure tested to 95kPa in accordance with BS EN 14401
- Leak tested in accordance with BS EN 14254

Choice of Volume with Wide Mouth Design

- 120ml or 90ml volume suitable for most applications
- Wide 53mm opening for ease of sample collection and management

Available both Sterile and Non-Sterile

 Sterilised by gamma radiation for Sterility Assurance Level (SAL) 10⁻⁶

CE Mark

 Compliant with the requirements of In Vitro Diagnostic Device Directive 98/79/EC

Excellent Clarity with — Easy to See Graduations

• For visibility of samples within the vial

Stackable

Optimize sample storage

Knurls on Cap and Base

Easy to handle with gloved hands

The term "leakproof" applies to Thermo Scientific Samco Clicktainer products, 90ml/53mm and 120ml/53mm, that meet the following criteria: a) a vacuum test that draws 95 kPa of pressure within a sealed chamber with no leakage; b) leak testing in accordance with BS EN 14254, which involves filling the container with a dye, placing for 2 minutes on a roller mixer followed by 1 hour inverted in a water bath with no leakage. Note: these tests, using other liquids, may not yield the same results. To ensure safe usage, customers are advised to test Thermo Scientific Samco Clicktainer Vials under conditions of their planned applications. Thermo Scientific Samco Clicktainer Vials are leakproof at ambient temperature and pressure when used with their Clicktainer closures.

ISO Compliant: Tested in accordance with the requirements of BS EN 14254 and BS EN 14401 standards



Label Design, Cap Colors and Formats

Clearly designed label for patient fields and required information

- Available with and without Label
- Temperature Strip Label
- Temperature Strip only
- One Label, Multiple Languages English, Spanish, French
- Required FDA and CE mark information on the label

Wide Range of Cap Colors and Packaging Formats

- Wide range of cap colors helps with sample differentiation and management
- Available as assembled product for general laboratory use, clinical applications and mice necropsy
- Available as separate components for manufacturing automation where suppliers add formalin solution, boric acid, sodium thiosulfate or other reagents for resale into histology, urinalysis, water testing, etc.
- Individual packaged product for sterile environments such as hospital operating rooms and cleanrooms
- laboratories and doctors offices







thermoscientific

Thermo Scientific Samco Clicktainer Vial



	Product Code	Description	Cap Color	Format	Inner Pack/Case Qty
120mL/53mm					
	15702197	SAMCO Clicktainer Vial, Tabbed Label, Sterile	Orange	Preassembled	75/300
	15722197	SAMCO Clicktainer Vial, Tabbed Label, Sterile	O White	Preassembled	75/300
	15732197	SAMCO Clicktainer Vial, No Label, Non-Sterile	O White	Separate Cap & Vial	300/300
	15762197	SAMCO Clicktainer Vial, Non-tabbed Label, Non-Sterile	Yellow	Preassembled	75/300
	15742197	SAMCO Clicktainer Midstream Urine Kit, No Funnel, Sterile	O White	Preassembled	1/100
	15712197	SAMCO Clicktainer Vial, No Label, Peel Pouch, Sterile	Orange	Preassembled	1/100
	15752197	SAMCO Clicktainer Vial, No Label, Peel Pouch, Sterile	O White	Preassembled	1/100
90mL/53mm					
	15782197	SAMCO Clicktainer Vial, Tabbed Label, Sterile	Orange	Preassembled	100/400
	15702207	SAMCO Clicktainer Vial, Tabbed Label, Sterile	O White	Preassembled	100/300
	15772197	SAMCO Clicktainer Vial, Tabbed Label, Green Cap, Sterile	Green	Preassembled	100/400
	15792197	SAMCO Clicktainer Vial, Tabbed Label, Sterile, Dispenser Pack	Orange	Preassembled	20/320
	15712207	SAMCO Clicktainer Vial, Temperature Strip Label, Sterile	O White	Preassembled	100/300
	15722207	SAMCO Clicktainer Vial, Temperature Strip Only, Sterile	O White	Preassembled	100/400
	15732207	SAMCO Clicktainer Vial, Non-Tabbed Label, Non-Sterile	O White	Preassembled	75/300
	15782207	SAMCO Clicktainer Vial, Non-tabbed Label Non-Sterile	Yellow	Preassembled	100/400
	15742207	SAMCO Clicktainer Vial, Temperature Strip Label, Non-Sterile	O White	Preassembled	300/300
	15752207	SAMCO Clicktainer Vial, No Label, Non-Sterile	O White	Separate Cap & Vial	300/300
	15762207	SAMCO Clicktainer Vial, No Label, Non-Sterile	O White	Preassembled	300/300
	15772207	SAMCO Clicktainer Vial, Temperature Strip Only, Non-Sterile	O White	Preassembled	100/400

Find out more at www.eu.fishersci.com/go/clicktainer

© 2017 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at fishersci.com/trademarks.

Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303 Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

