

Thermo Scientific Value Lab Refrigerators, Freezers and Combination Units

An economical solution for safe, sample storage.



3.2 cu. ft. lab refrigerator

Thermo Scientific™ Value laboratory refrigerators, freezers and combinations units are particularly suited for routine laboratory applications with sizes ranging from small benchtop to full-height freestanding. Our features are not available in domestic models and meet the demanding requirements of the laboratory environment.

Available in North America and select Latin America countries, consult your sales representative.

Alert: These units are not intended for critical lab samples. Storage of high-value, life science materials or clinical biologicals may require alarm, monitoring,

back-up and temperature performance features. The need to add these features should be determined by your industrial hygienist or safety officer. These units are not suitable for corrosive environments.

- Factory tested prior to shipment
- NRTL listed
- Environmentally-friendly, non CFC refrigerants
- Adjustable temperature control
- Easy-to-clean, durable finish
- Manual and auto defrost combination units
- Manual defrost refrigerators and freezers

Warranty: 13 months parts and labor

Certifications: CSA and or UL



10 cu. ft. combination lab refrigerator and freezer



20.6 cu. ft. lab freezer



24.6 cu. ft. combination lab refrigerator and freezer

Thermo Scientific Value Lab Refrigerators (Temperature Range 1° to 12°C)

Cat. No.	Capacity cu. ft.	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
14-512-520	1.8	manual	115V / 60Hz	1 fixed, 2 door	yes	13 x 16.5 x 15.75**	21 x 19.75 x 20.75	55
13-990-530*	5.5	manual	115V / 60Hz	3 fixed, 2 door	yes	15 x 20 x 26.5	23.5 x 23.75 x 33.75	110
13-989-19	5.6	manual	115V / 60Hz	4 adjustable	yes	16 x 20 x 27.5***	25.13 x 23.65 x 33.5	102
14-512-526*	6.1	manual	115V / 60Hz	2 fixed, 3 door	no	15 x 21 x 26.5†	24 x 24 x 34	120
11-651-121	20.3	manual	115V / 60Hz	5 fixed	yes	20 x 27 x 57††	32 x 30.75 x 70	240

* Features internal digital temperature display † 8" (for bottom 5"), 15" (for top 21.5")

** 8" (for bottom 8"), 13" (for top 8") †† 17" (for bottom 20"), 20" (for top 40")

*** 10" (for bottom 8"), 16" (for top 19.5")

Thermo Scientific Value Lab Freezers (Temperature Range -12° to -20°C)

Cat. No.	Capacity cu. ft.	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
14-512-521	1.6	manual	115V / 60Hz	1	yes	13.5 x 15.5 x 14*	21 x 19.75 x 20.75	55
14-512-523**	3.2	manual	120V / 60Hz	3 drawers	yes	16.5 x 17.75 x 24.25	23.5 x 23.5 x 33.5	102
3752-DB	5	manual	115V / 60Hz	2 fixed, 3 door	no	15 x 18 x 28.5	23.5 x 21 x 33.5	84
11-675-951	20.6	manual	115V / 60Hz	4 fixed	yes	20 x 27 x 57	31 x 32 x 71	200

* 7" (for bottom 8"), 13.5" (for top 6")

** -15 to -35°C, suitable for -25°C applications; model contains a high-temperature alarm

Thermo Scientific Value Combination Lab Refrigerator/Freezers (Temperature Range: 1° to 12°C [ref], -12° to -20°C [frz])

Cat. No.	Capacity cu. ft.	Defrost	Voltage	Shelves Interior/Door	Drawers	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
22-299-773	5.6	manual	120V / 60Hz	3 / 3	–	18 x 15 x 21 (ref.) 11 x 15 x 5 (frz.)	25.5 x 21.25 x 33.5	78
14-512-524	8.8	automatic	115V / 60Hz	2 / 2	1	13 x 19.25 x 33 (ref.) 8.75 x 18.5 x 14.5 (frz.)	23.5 x 23.25 x 60.75	160
11-651-120	10	manual	120V / 60Hz	Ref: 3 / 4, Frz: 1	–	13.5 x 19.5 x 38.5 (ref.) 15 x 17.5 x 11.5 (frz.)	26.5 x 23.5 x 60	150
14-512-522†	12	manual	115V / 60Hz	Ref: 4	Frz: 3	13.5 x 19 x 35 (ref.) 13.5 x 17 x 28 (frz.)	23.75 x 23.75 x 79	251
18LC-16WW-TS	18.2	automatic	120V / 60Hz	Ref: 2 full Frz: 1 half, 2 full	Ref: 3	18.75 x 22.5 x 39.5 (ref) 14.5 x 23.75 x 13.5 (frz)	32.4 x 29.6 x 66.1	205
13-991-90**	24.6	automatic	120V / 60Hz	Ref: 4 adj. half 1 full / 6 door Frz: 1 / 3	3	20 x 32 x 37 (ref.) 17 x 31 x 18.5 (frz.)	36 x 35.75 x 68	242

* Key locks are not available on combination units

** Features internal digital temperature display

† Dual compressor system, automatic refrigerator and manual defrost freezer



© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific